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
口蹄疫病毒样颗粒组装机制及免疫特性的研究  
Research on the assembly mechanism and immune characteristics of foot-and-mouth disease virus-like particles

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
研究背景  
Research Background

Science-driven solutions®

**病毒样颗粒** (Virus-like particles, VLPs) 是一种由病毒结构蛋白组成但不含病毒核酸、不能自主复制的颗粒。 Virus-like particles (VLPs) are particles composed of viral structural proteins but without viral nucleic acids and unable to replicate autonomously.

与天然病毒粒子更接近  
Closer to natural virus particles

Monolayer

1

3

VLPs具有比可溶性蛋白更强的免疫原性  
VLPs have stronger immunogenicity than soluble proteins

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




### 技术优势 Technical advantage

- ◆ 无需高致病性生产资质
- ◆ No need for high pathogenicity production qualification
- ◆ P1-2A聚合体原核表达困难，表达量低，且P1蛋白后期加工存在技术瓶颈
- ◆ Prokaryotic expression of P1-2A polymer is difficult, with low expression levels, and there are technical bottlenecks in the later processing of P1 protein
- ◆ 单个VP1蛋白体外以包涵体形式表达
- ◆ Expression of single VP1 protein in vitro in the form of inclusion bodies

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"Foot-and-mouth disease virus-like particle vaccine: a leading generation of prevention and control tool"

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The foot-and-mouth disease virus-like particle vaccine is the world's first virus-like particle vaccine self-assembled from three structural proteins in a cell-free system, and has obtained two Class I new veterinary drug certificates.

猪口蹄疫O型病毒样颗粒疫苗一类证书  
Class I Certificate for Porcine Foot-and-Mouth Disease Type O Virus-Like Particle Vaccine

牛口蹄疫O型病毒样颗粒疫苗一类证书  
Class I Certificate of Foot-and-Mouth Disease O-type Virus-like Particle Vaccine for Cattle


猪口蹄疫O型、A型二价VLPs疫苗、牛口蹄疫OA二价VLPs疫苗。  
Bivalent VLPs vaccine against porcine foot-and-mouth disease type O and type A, and bivalent VLPs vaccine against bovine foot-and-mouth disease type O and type A.

PART 2

口蹄疫VLPs的组装  
Assembly of Foot-and-Mouth Disease VLPs

Science-driven solutions®

01

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
Foot-and-mouth disease

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Quantitative analysis of toxic  
particles

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Figure 5. Flowchart of FMD VLPs isolation

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Figure 1. Results of sucrose density gradient centrifugation of VLPs

VLPs的出峰时间稍晚与FMDV，由于人工组装的特点，VLPs的峰形比较舒缓。The peak time of VLPs is slightly later than that of FMDV. Due to the characteristics of artificial assembly, the peak shape of VLPs is relatively smooth.

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Figure 2. Identification results of VLPs collected by sucrose density gradient centrifugation

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The collection peak of VLPs showed three bands of foot-and-mouth disease virus, and DLS and electron microscopy detected particles with a size of about 30 nm

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Figure 6 Analysis results of VLPs and monomeric protein antibody levels

VLPs诱发的免疫反应更持久，可诱发高水平特异性抗体和中和抗体；单体诱发的免疫反应持续期短，且诱发中和抗体的能力不足。  
The immune response induced by VLPs is more durable, capable of eliciting high levels of specific and neutralizing antibodies; the immune response induced by monomers is short-lived and lacks the ability to induce neutralizing antibodies.

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
Figure 1. Assembly model of foot-and-mouth disease virus-like particles

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Figure 8 Effect of assembly solution pH on VLPs assembly

A. SDS-PAGE analysis results; B. DLS analysis results; C. VLPs content detection results

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口蹄疫病毒pH值耐受范围较窄，而病毒样颗粒的组装试验显示：组装液在pH8.0的条件下组装率明显高于pH8.5  
The pH tolerance range of foot-and-mouth disease virus is relatively narrow, and the assembly experiment of virus-like particles shows that the assembly rate of the assembly solution at pH 8.0 is significantly higher than that at pH 8.5

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Figure 10 Effect of dispersant on the assembly of O-type VLPs

A. SDS-PAGE results; B. DLS results; C. VLPs content detection results; D. Effect of E reagent concentration on VLPs assembly rate  
A. SDS-PAGE结果; B. DLS结果; C. VLPs含量检测结果; D. E试剂浓度对VLPs组装率的影响

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Adding dispersing agents to purified proteins can significantly increase the assembly rate of O-type VLPs

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Figure 11 Effect of dispersant on the assembly of type A VLPs

A. Turbidity comparison effect diagram; B. DLS results; C. VLPs content detection results

在纯化蛋白中加入分散试剂亦可明显提高A型VLPs的组装率，且A型口蹄疫病毒结构蛋白的聚集明显比O型严重  
Adding dispersing agents to the purified protein can also significantly increase the assembly rate of type A VLPs, and the aggregation of type A foot-and-mouth disease virus structural proteins is significantly more severe than that of type O

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Figure 13 Analysis results of the effect of F reagent on the assembly of O-type VLPs

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F reagent can stabilize the structure of VLPs, and when used together with N reagent, it can increase the assembly rate of VLPs

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- 在**组装后的VLPs**中添加不同浓度的B试剂，放置16h后检测，VLPs含量0.05%>0.03%>0
- Add different concentrations of B reagent to the assembled VLPs, and detect after 16 hours. The VLPs content is 0.05%>0.03%>0

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

Figure 14 Analysis results of the effect of B reagent on the assembly of O-type VLPs

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The results showed that the addition of B reagent to VLPs significantly increased the content of VLPs

## Stabilization of Human Papillomavirus Virus-Like Particles by Non-Ionic Surfactants

LI SHI, GAUTAM SANYAL, ALEX NI, ZHENG LUO, SARAH DOSHNA, BEI WANG, TAMMY L. GRAHAM, NING WANG, DAVID B. VOLKIN

Biologics and Vaccines, Pharmaceutical Research and Development,  
Merck Research Laboratories, West Point, Pennsylvania 19486

Received 7 October 2004; revised 9 February 2005; accepted 28 February 2005

Published online in Wiley InterScience (www.interscience.wiley.com). DOI 10.1002/jps.20377

**ABSTRACT:** Human papillomavirus (HPV) virus-like-particles (VLPs) produced by recombinant expression systems are promising vaccine candidates for prevention of cervical cancers as well as genital warts. At high protein concentrations, HPV VLPs, comprised of the viral capsid protein L1 and expressed and purified from yeast, are protected against detectable aggregation during preparation and storage by high concentrations of NaCl. At low protein concentrations, however, high salt concentration alone does not fully protect HPV VLPs from aggregation. Moreover, the analytical analysis of HPV VLPs proved to be a challenge due to surface adsorption of HPV VLPs to storage containers and cuvettes. The introduction of non-ionic surfactants into HPV VLP aqueous solutions provides significantly enhanced stabilization of HPV VLPs against aggregation upon exposure to low salt and protein concentration, as well as protection

- VLPs在组装或放置过程中会出现聚集，聚集会使组装率下降，且会使VLPs不稳定
- VLPs may aggregate during assembly or placement, which can reduce the assembly rate and destabilize the VLPs
- 及免疫原性下降或丧失。 and decreased or lost immunogenicity.

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Table 1: Statistical results of immune challenge protection of porcine foot-and-mouth disease virus type O VLPs vaccine

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The results of the antiviral protection test showed that aggregation reduces the immunogenicity of VLPs, while inhibiting aggregation can enhance the immunogenicity of VLPs, but excessive inhibition can destroy the immunogenicity.

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




The results of the freeze test showed that there was no significant difference in the concentration of VLPs stored at -20°C and -70°C for 5 days compared to those stored at 4°C. Additionally, the peak overlap of the sucrose density gradient was good, indicating that VLPs have a certain degree of tolerance to freezing.

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
The mechanical shear results showed that there was no significant difference in VLPs concentration between VLPs sheared at 3000 rpm/min and 300 rpm/min for 15 minutes and VLPs stored at 4°C. The peak overlap of the sucrose density gradient was good, indicating that VLPs had certain mechanical resistance.

PART 3

口蹄疫VLPs的免疫特点  
Immunological characteristics of foot-and-mouth disease VLPs

Science-driven solutions®

01

口蹄疫VLPs疫苗免疫力的产生  
Immunity generation of  
foot-and-mouth disease  
VLPs vaccine

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Table 1 Statistics of antibody and challenge incidence after immunization with foot-and-mouth disease VLPs vaccine at different times

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After emergency immunization, the antibody positive rates and protection rates at 3 days, 7 days, and 14 days were 0%, 0%, 60%, and 0%, 40%, and 100%, respectively.

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

The 3ABC antibody test results after virus challenge showed that the 3ABC antibodies of the non-infected pigs were all negative, indicating that the virus was completely neutralized and could not replicate.

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Following an emergency immunization challenge, the body temperature of animals immunized for 14 days remained normal throughout the entire challenge period. Although the 7-day immunization group only achieved 40% protection, their body temperature was significantly lower than that of the control group and the 3-day immunization group.

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
Figure 13 Changes in FMDV content in serum and tissues after emergency immunization

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After emergency immune challenge, the serum and tissues of animals immunized for 14 days were all negative for FMDV; although only 40% of animals in the 7-day immunization group were protected, the virus content in serum and tissues was significantly lower than that in the control group and the 3-day immunization group.


02

口蹄疫VLPs疫苗免疫特点及持续期  
Immunization characteristics  
and duration of foot-and-mouth  
disease VLPs vaccine

A

B


图1 猪口蹄疫O型、A型二价VLPs疫苗与猪全病毒灭活疫苗抗体水平比对结果 (猪)  
Figure 1 Comparison of antibody levels between bivalent VLPs vaccine against foot-and-mouth disease type O and type A and inactivated whole virus vaccine in pigs

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The results of antibody level detection after immunizing pigs with VLPs vaccine and whole virus inactivated vaccine showed that the production period of antibodies induced by VLPs vaccine was significantly shorter than that of inactivated vaccine, with a specific antibody level of 1:64 achieved within 14 days; the neutralizing antibody level was not lower than that of inactivated vaccine.

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Figure 13 Comparison of the change trend of O-type antibody between VLPs vaccine and high-end inactivated vaccine

VLPs疫苗诱发抗体的产生速度要明显快于灭活疫苗，14天70%的抗体转阳，且单次免疫在5个月内抗体水平均不低于灭活疫苗。The production speed of antibodies induced by VLPs vaccine is significantly faster than that of inactivated vaccine, with 70% of antibodies turning positive within 14 days, And the antibody level of single immunization is not lower than that of inactivated vaccine within 5 months.

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Table 5 Statistical table of antibody level and challenge of swine foot-and-mouth disease type O VLPs vaccine

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The antibody persistence period of the porcine foot-and-mouth disease virus-like particle vaccine after a single immunization is 6 months, and the protection rate after 6 months of immunization is greater than or equal to 80%.

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Bivalent virus-like particle vaccine for foot-and-mouth disease of cattle, types O and A: Each dose contains 13.9 PD<sub>50</sub> for type O and 15.5 PD<sub>50</sub>

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
猪口蹄疫O型、A型二价病毒样颗粒疫苗: O型每头份含14.3 PD<sub>50</sub>; A型每头份含14.3 PD<sub>50</sub>  
Bivalent virus-like particle vaccine against foot-and-mouth disease in pigs, types O and A: Each dose contains 14.3 PD<sub>50</sub> for type O and 14.3 PD<sub>50</sub> for type A

03

口蹄疫VLPs疫苗稳定性  
Stability of foot-and-mouth  
disease VLPs vaccine


图8 放置9个月VLPs疫苗的抗体检测结果  
Figure 8 Antibody detection results of VLPs vaccine placed for 9 months

免疫试验结果显示：放置9个月的口蹄疫病毒样颗粒疫苗免疫60天的架子猪，  
 单次免疫30天后抗体水平合格率100%。 The results of the immune test showed that the 60-day-old pigs immunized  
 with the foot-and-mouth disease virus-like particle vaccine placed for 9 months,  
 The qualified rate of antibody level was 100% after 30 days of single immunization.

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Table 11 Statistical results of immune challenge protection of foot-and-mouth disease type O VLPs vaccine

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The results of the virus attack protection test showed that the foot-and-mouth disease virus like particle vaccine placed for 9 months can still protect target animals against the attack of the O-type virulent strain, and the protective efficacy still meets the standard of PD<sub>50</sub> greater than 6 in China.

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
Efficacy test of vaccine after 18 months of storage

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After being stored for 18 months, the efficacy of both batches of vaccines decreased, but the PD<sub>50</sub> of both batches of vaccines remained greater than 6.

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Four types of foot-and-mouth disease VLPs vaccines have been launched on the market. VLPs vaccines are a new generation of vaccines with broad antigenicity, safety, high efficiency, and precise differentiation between infected and immune animals. They are specifically designed for the purification of foot-and-mouth disease.

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We kindly request that all leaders and experts provide their critiques and corrections, and continue to monitor the FMD VLPs vaccine

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